

Corncob Twist

Teklad 7076



NEW! Sustainable enrichment bedding

PRODUCT DESCRIPTION

Teklad 7076 Corncob Twists is the newest bedding product that combines our two popular products - corncob bedding with 1/2 inch pieces of Diamond Twists.

Introducing the ultimate bedding and enrichment solution - the perfect combination of economical bedding and enrichment! Our all-natural and sustainable corn cob bedding, with its superior absorbency and ammonia control, are made from the woody-ring portion of the corncob without any other components or additives. During production, whole corncobs are ground, dried, sifted and aspirated to make very uniform particles that are essentially dust free. The entire process is computer controlled and monitored and includes regular quality control measurements for moisture and particle size distribution.

For added animal welfare benefit and enrichment, we added our Diamond Twist enrichment product. These 1/2" pieces of thick rolled tissue paper offer low dust, ensuring a clean and healthy environment. This product is shreddable and does not wind around rodents legs or pups. The paper material will not clog drains or disposals.

PRODUCT FEATURES AND BENEFITS

- Eliminates the need for added enrichment.
- Encourages natural foraging behavior in rodents
- Improves air flow and is highly absorbent
- Available in irradiated and autoclavable packaging
- Available in 1/8" and 1/4" sized corncob

PACKAGING

Corncob Twist is packaged in:

- 28 lb irradiated case (4 x 7 lb bags)
- 25 lb irradiated bag
- 40 lb autoclavable perforated bag

Find out more
Contact your Inotiv representative
or visit us at inotivco.com