

TD.130856

TAM Diet (250, 2016, Red)

Formula

	g/Kg
2016, Teklad Global 16% Protein Rodent Diet	949.75
Sucrose	49.75
Tamoxifen USP	0.25
Red Food Color	0.25

Footnote

Formulated for 250 mg tamoxifen/kg diet. This will provide ~40 mg tamoxifen per kg body weight per day assuming 20-25 g body weight and 3-4 g intake. Dose consistent with Welle, et. al. AJPEM (2009) 297(4):E935-40. 2016 base provides a minimal phytoestrogen background. Contains red food dye as a visual indicator. Tamoxifen is premixed with sucrose however feed aversion is still common with tamoxifen containing diets. See <http://www.inotivco.com/tamoxifen-diets> for more information.

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	15.4	19.7
Carbohydrate ²	55.1	70.6
Fat	3.4	9.7

Kcal/g 3.1

¹ Values are calculated from ingredient analysis or manufacturer data

² Estimated digestible carbohydrate

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Standard Diet Base
- Teklad Global 2016
- Tamoxifen Inducible Cre
- Color Coded Red

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com